

ADAPTERY DO NAWILŻACZA BUTLOWEGO



RespiFlo 34x1,5



RespiFlo 16x1

NAZWA

gwint do
dozownika

gwint do
nawilżacza

Adapter RespiFlo 14x1 wew. (stary FARUM)	14x1 wew.	9/16	NBZ 01 08 00
Adapter RespiFlo 34x1,5 zew. (DRAGER)	34x1,5 zew.	9/16	NBZ 01 09 00
Adapter RespiFlo G3/8 wew. (francuskie)	G3/8 wew.	9/16	NBZ 01 10 00
Adapter RespiFlo 16x1 zew. (nowy FARUM)	16x1 zew.	9/16	NBZ 01 12 00
Adapter RespiFlo 12x1,25 wew. (włoskie)	12x1,25 wew.	9/16	NBZ 01 13 00
Adapter RespiFlo G1/4 wew. (korgiel)	G1/4 wew.	9/16	NBZ 01 14 00
Adapter RespiFlo 9/16 wew.	9/16 wew.	9/16	NBZ 01 15 00

ADAPTER BUTLOWY FARUM



NAZWA

Adapter butlowy do dozownika FARUM



INDEKS

NBZ 01 05 00

TULEJA



Służy do bezpośredniego podłączenia węża do dozownika.

INDEKS

TN 01 02 10